Operations Guide

MODEL 502
MODEL 503
MODEL 504
MODEL 505
MODEL 506
MODEL 507

QHL75[®]

Curing Light

POLYMERIZATION UNIT

CE

TECHNICAL DETAILS

English

AC INPUT RATING: MODEL 506: 120V, 50/60Hz, 2.0A
 MODEL 502, 503, 504, 507: 230V, 50/60Hz, 1.0A
 MODEL 505: 100V, 50/60Hz, 2.0A

FUSE RATING: MODEL 506: T2.5A, 250V
 MODEL 502, 503, 504, 507: T1.6A, 250V
 MODEL 505: T2.5A, 250V

LAMP RATING: 75W, 11V

Protection from Electric Shock: Type BF Applied Part
 IEC Equipment Class I

Protection against ingress of water: Ordinary Equipment

Operation Mode: Continuous operation with intermittent loading:
 (1 min./5 min. = 1 min. ON / 5 min. OFF)

Operating Environment: **NOTE: Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with nitrous oxide.**

Temperature: 10 to 25 deg. Celsius
 Relative Humidity: 30% to 75% (non-condensing)
 Atmospheric Pressure: 700 to 1060 hPa

Transport and Storage Conditions:

Temperature: -40 to +70 deg. Celsius
 Relative Humidity: 10% to 100% (non-condensing)
 Atmospheric Pressure: 500 to 1060 hPa

Dimensions:

Height: (base):	3.62 inches	(9.18 cm)	Hand-piece cord		
(with hand-piece):	9.12 inches	(23.15 cm)	length (approx.):	80.0 inches	(203.0 cm)
Width:	6.00 inches	(15.22 cm)	Power cord length:	120 inches	(305.0 cm)
Length:	6.13 inches	(15.55 cm)	Weight:	5.78 pounds	(2.63 kg)

NOTE: - Models 502, 503, 504 and 506 meet all of the following Safety Standards.
- Model 505 is pending approval certifications.

Complies to following Directives/Standards:

CE - 93/42/EEC	- MDD Directive
CE - 89/336/EEC	- EMC Directive
IEC 601-1	- 1988 - Medical Electrical Equipment (General Requirements for safety)
IEC 601-1 Am.1	- 1991 - Amendment 1 to IEC 601-1
IEC 601-1 Am.2	- 1995 - Amendment 2 to IEC 601-1
EN 60601-1	- 1993 - Medical Electrical Equipment (General Requirements for safety)
EN 60601-1-2	- 1993 - IEC 601-1-2: Medical electrical equipment Part 1: General requirements for safety 2: Collateral standard: Electromagnetic compatibility - Requirements and tests

NOTE: A QHL75® Curing Light Repair Manual is available upon request from DENTSPLY.

GLOSSARY OF SYMBOLS

Attention, Consult Accompanying Documents

Alternating Current

Direct Current and Alternating Current

Protective Earth

OFF (Power: Disconnection from the Mains)

ON (Power: Connection to the Mains)

Type BF Applied Part

INTRODUCTION

The QHL75® Curing Light offers the combination of features, performance and reliability for the polymerization of all photosensitive dental materials (in the 400 to 500nm visible light range.) This versatile dental curing light system offers:

- High Light Output Power (minimum 450 mW/cm²).
- Cycle ON/OFF trigger switch with cycle interrupt capability.
- Choice of 5 standard DENTSPLY Light Probes for all resin polymerization needs.
- Optional Wall Mounting Kit available.

The QHL75® curing Light will assure your complete satisfaction for all resin polymerization needs.

INSTALLATION

Note: - Upon receiving the QHL75® Curing Light, check the packaging and parts for any possible damage that may have occurred in transit. If damage is apparent, please contact your dealer immediately.

- Verify the line plug of the power cord is properly rated, that all connections have been performed to safe practices and all local regulation requirements (always use the light in conjunction with a properly grounded outlet).
- Affix the eye protection shield to the probe chuck on the nosecone of the hand-piece or use suitable protective eye wear.

1. Insert the desired fiberoptic probe into the nozzle of the handpiece until it snap-locks in place.
2. Insert the holster in mounting hole on top of the base/power unit (forks positioned forward) and place the handpiece in holster with the fiberoptic probe/nosecone facing the rear of unit.
3. Plug the detachable power cord into receptacle on the side of the base/power unit. The power cord may now be plugged into the AC supply.
4. Refer to the CLEANING section for specific asepsis/disinfecting instructions.

OPERATION

The main power switch is located on the side of the base/power unit. The switch illuminates (green) when pressed to the "on" - "I" position (Models 502 - 506) indicating the unit is receiving power and is ready for curing. The main power is now controlled by the activation (trigger) switch on the handpiece.

Momentarily depress the trigger switch on the front of the handle to activate the lamp and begin the exposure cycle. The unit will beep every 10 seconds during the exposure.

The QHL75® Curing Light will automatically shut OFF after 1 minute of exposure time. The unit can be activated again by momentarily depressing the on/off trigger switch on the front of the handle.

The exposure cycle can be interrupted at any time by depressing the trigger switch once again.

The high speed cooling fan will start at the beginning of each exposure. After the final exposure, the fan will continue to run for approximately five to six minutes before automatically shutting off. **Note: Do not turn off main power or unplug the unit until the fan has finished this cycle. The QHL75® Curing Light is protected from over-heating by thermal safety devices. Should the unit overheat, it will automatically turn off until normal operating temperatures are reached.**

When the curing operation is complete, return the handpiece to the handle holster. During prolonged periods of non-use, the unit may be turned off by depressing the main power switch to the "off" - "O" position, switch illumination will also turn off (Models 502-506).

CURING:

- For optimal curing, always refer to our composite manufacturer's exposure time recommendations.
- **Fiberoptic Light Probes:** The QHL75® Curing Light comes with one 8mm, 60° probe. Additional probes may be ordered separately (Figure 5). All light probes rotate 360° within the probe chuck.
- Check and measure the light output intensity daily using a light intensity radiometer to verify lamp, filter, probe performance and to ensure proper curing of resin materials. The **CureRite™** by DENTSPLY is an excellent radiometer.

Attention: Frequently verify the light output intensity (using the standard 8mm fiberoptic probe) on a CureRite™ radiometer to ensure sufficient curing power.

CONTRAINDICATIONS

1. Persons who have a medical history indicating photobiologic reactions (including persons with solar urticaria or erythropoietic protoporphyria) or who are using photosensitizing drugs (including 8-methoxypsoralen or tetracycline) should not be exposed to light from the lamp.
2. In cases where a patient or user of the QHL75® Curing Light has been fitted with an implanted heart pacemaker and has been cautioned against the use of small electrical appliances (such as electric shavers, hair dryers, etc.) we recommend that the QHL75® Curing Light not be used.

WARNING

1. When curing the material, position the light over the polymerizable material and avoid extraneous contact with soft tissue whenever possible. Prolonged continuous exposure to soft tissue can result in injury or irritation.
2. Allow the QHL75® Curing Light maximum time to cool prior to removing the nosecone, exposing potentially hot metal lamp reflector cups.

SAFETY PRECAUTIONS

1. Never activate lamp when the fiberoptic probe has been removed from handpiece.
2. Do **NOT** look directly at the light emitted from the tip of the instrument. As with any bright light source, staring directly at the source may cause temporary effects such as after-image.
3. A filter is located within the handpiece to allow transmission of only useful curing irradiation. The filter should be examined before use for any cracks or defects.

A CRACKED OR MISSING FILTER MUST BE REPLACED BEFORE USE TO PREVENT THE EMISSION OF HARMFUL IRRADIATION.

4. Persons with known retinal disease should be careful and should employ suitable protective measures after consultation with their ophthalmologist.
5. Persons who have had surgery for eye cataracts are especially sensitive to light and generally should not use the QHL75® Curing Light unless they employ suitable protective measures which should include the use of spectacles which filter out blue-violet and ultraviolet light.
6. When a lamp has burned out wait until the lamp is cool before attempting to replace it.
7. Fiberoptic tips must be maintained in a clean condition. Contaminating the tip will reduce the effectiveness of the light transmitted.
8. Other lamps may appear similar but perform differently. Replace only with genuine QHL75® Curing Light replacement bulbs.
9. DO NOT IMMERSE UNIT IN WATER. See "Cleaning Instructions".

CLEANING

The fiberoptic light probe may be autoclaved or dry heat sterilized (see manufacturer's instructions). Any cured resin materials on the tip of the light probes should be removed (pried off) with a razor blade or sharp knife to ensure proper curing of composite materials.

NOTE:

Always protect open ends of the light probes from abrasive materials. Should polished ends become marred, light transmission and curing will be diminished.

Disconnect the power cord from the base unit prior to cleaning.

Do not use a glutaraldehyde based solution to disinfect your unit. This could affect the mechanical operation and damage the plastic housing.

Clean and disinfect the surfaces of the unit housings by applying an approved non-immersion type disinfectant solution* carefully following the instructions provided by the disinfectant solution manufacturer.

*NOTE: Water based disinfectant solutions are preferred. Alcohol-based disinfectant solutions may be harmful and may discolor plastic materials.

Perform the disinfectant procedure only after the unit has cooled to room temperature. Wipe the unit with a damp cloth to remove any residual deposits of the disinfectant solution.

Use of alcohol-based disinfectant solutions will void the warranty.

Remember to reconnect the power cord to the unit when cleaning process is completed.

LAMP REPLACEMENT

The lamp used in the QHL75® Curing Light is a 75 watt, 11.0 volt quartz halogen type, which has an estimated life rating of 100 hours under normal usage and with nominal input voltages (increased or decreased AC input voltage may reduce life). Bulb life may be shortened by a sharp blow during use or cool down period. As the light intensity of all bulbs decreases during usage, we recommend to check the light output frequently using a **CureRite™**.

Replacement is recommended when lamp transmission is below the minimum specified by composite manufacturers or less than a 300mW/cm² radiometer reading.

Disconnect the power cord from the unit and be sure the lamp has cooled before replacing.

Proceed by turning the nosecone to align the line with the open circle on the body of the handpiece. Pull the nose cone off to expose the light bulb. Gently wiggle the lamp out of the socket and replace with a new lamp. A bulb extraction tool is provided should you have difficulty removing the bulb.

Good practice is to inspect the optical filter within the nose cone. Remove any contaminants using a cotton swab. A cracked, peeling, discolored or defective filter should be replaced.

When inserting the new lamp, it is recommended that you use a tissue to hold the lamp and prevent skin acids from etching the quartz envelope which will lead to premature lamp failure. Replace the nose cone and turn so the line is opposite the closed circle on the handpiece

FUSE REPLACEMENT

NOTE:

Always disconnect the power cord from the unit before replacing the fuses (use fuses with the same ratings).

1. The fuses are located in the fuse compartment of the AC receptacle, at the side of the base unit.
2. With a small screwdriver, pull out the fuse compartment and let the door fall open (on its pivot) to expose the fuses.
3. Pull out the fuses, visually inspect the fuses and replace if necessary with fuses of the same ratings.
4. Lift the fuse door open until it covers the fuses and push the fuse compartment completely in the AC receptacle until it locks in place.

Remember to reconnect the power cord.

WARRANTY

This product is designed for use in a dental office and this warranty is not applicable to other uses. This product is warranted against defective material and workmanship for a period of twenty-four (24) months (with an additional thirty-six (36) months warranty on the handpiece cord) from the date of purchase by a dentist. In the event of such a defect, DENTSPLY will repair or replace the product or necessary parts therein, at its discretion, and such repair or replacement shall be the sole remedy of this warranty. This warranty extends only to the original purchase and is subject to these conditions.

- 1/ Any servicing of this product must be performed by qualified DENTSPLY Dealer service personnel. The product must be returned in its original packaging or charges may apply.
- 2/ The unit must not be subjected to abuse or improper installation or application, nor installation or repair service by other than trained DENTSPLY service personnel.
- 3/ Any light bulbs or fiberoptic light probes are excluded from this warranty.

THERE ARE NO WARRANTIES, EXPRESS OR LIMITED, WHICH EXTEND BEYOND THIS DESCRIPTION. DENTSPLY neither assumes, nor authorizes any person to assume for it, any other liability in connection with the sale and use of this product.

DAMAGES ARE LIMITED STRICTLY TO REPLACEMENT OF THE PRODUCT. DENTSPLY EXPRESSLY DISCLAIMS LIABILITY FOR INCIDENTAL AND CONSEQUENTIAL DAMAGES RESULTING FROM IT'S USE.

Claims covered by this warranty will be honored when this warranty is presented through your DENTSPLY Dealer within one (1) month from discovery of defect.

ACCESSORIES AND SPARE PARTS

Description:	European Reorder No.	Reorder No.
QHL 75® Replacement Bulb 75W/11,0V	603.17.507	644657
Wall Mounting Kit	603.21.200	644658
Curing probe, 13mm straight (bonding large facial-surfaces)	603.17.511	644710
Curing probe, 8mm 90° angle (posterior occlusal curing)	603.17.512	644711
Curing probe, 8mm 60° angle (general curing)	603.17.513	644712
Curing probe, 3mm 90° angle (incremental curing of posterior restorations)	603.17.514	644713
Curing probe, 13mm 80° angle (posterior occlusal curing)	603.17.515	644714

WORLDWIDE SERVICE CENTERS

DENTSPLY Professional
Attn.: Product Service
1301 Smile Way
P.O. Box 7807
York, PA 17404-0807, USA
1-800-989-8826

EU-Authorized Representative
DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
Germany
Tel. (49) (07531) 583-0

DENTSPLY Limited
Hamm Moor Lane
Addlestone
Weybridge, Surrey KT15 2SE
Tel. (44) (01932) 853422

DENTSPLY DeTrey
ZA du pas du Lac
17, rue Faraday
F-78180 Montigny le Bretonneux, France
Tel. (33) (01) 30147777

DENTSPLY Italia s.r.l.
Via A. Cavaglien, 26
I-00173 Roma, Italy
Tel. (39) (06) 7233626

Swiss Representative
DENTSPLY DeTrey Sàrl
Baar Office
Oberdorstr. 11
6342 Baar
Switzerland
Tel. (41) (041) 7662066

DENTSPLY Australia Pty Ltd.
204-206 Gipps Street
Abbotsford, 3067, Australia

DENTSPLY Japan K.K.
Tsunashima No. 2 Building
20-12 Yushima 3-chome
Bunkyo-ku, Tokyo 113, Japan

Hong Kong Division of DENTSPLY International Inc.
23/F Gee Ghang Hong Center
65 Wong Chuk Hang Road
Aberdeen, Hong Kong

DENTSPLY Mexico, S.A. de C.V.
Calzada Vallejo No. 846
Col. Industrial Vallejo
02300 Mexico, D.F.

DENTSPLY Industria e Comercio Ltda.
Rua Alice Herve 86
Bengen
25665-010 Petropolis, R.J., Brazil

DENTSPLY Caulk
Mfg by DENTSPLY Caulk
DENTSPLY International
Milford, DE 19963-0359

Manufactured by
DENTSPLY Caulk
Milford, DE 19963-0359
USA

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DENTSPLY Messico, S.A. de C.V.