The Palodent® System

Sectional Matrix System

DENSPLY
CAULK

DIRECTIONS FOR USE – ENGLISH

Caution: For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

The Palodent® System employs a spring steel BiTine® ring and sectional matrices for placement of restorations in the posterior region. Sectional matrices are available in three sizes suitable for use with all posterior restorative materials. BiTine® rings are available in round and oval or elongated shapes suitable for use in the posterior region either singly or in tandem. The advantages of the sectional matrices and rings include natural contours for better control of contact areas and embrasures, ease of placement, better visualization of the operative field, and more comfort for the practitioner and patient. Eliminating the leverage of a conventional matrix retainer allows the practitioner to use less tension, avoiding excessive compression of tooth and easier formation of anatomical proximal contact areas.

1.1 Delivery forms

- Sectional matrices are available in 3 sizes: Mini, Standard, and Plus. These matrices are suitable for use with all posterior restorative materials.
- BiTine® rings are available in 2 sizes, Standard and BiTine® II oval suitable for use in the posterior region either singly or in tandem.

1.2 Composition

- Sectional Matrices: Polished stainless steel
- BiTine® rings: Nickel plated high carbon spring steel
- Placement forceps: 420 Stainless steel

1.3 Indication

The Palodent® System is indicated whenever a restoration will involve interproximal surfaces to facilitate the contouring of a direct placement of restorative materials.

1.4 Contraindications

The sectional matrices are intended for single use only. Do not attempt to sterilize and/or re-use.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. The Palodent® System consists of small, sharp parts. As with all intraoral instruments, care must be exercised to ensure stable handling and placement. Placement of rubber dam or other suitable isolation is advised prior to operative procedures. Wear suitable protective eyewear and gloves.
Aspiration: If aspiration of matrix or BiTine® ring occurs seek immediate medical attention.
Ingestion: If accidental swallowing occurs, encourage patient to seek medical assistance to locate swallowed part. If nausea or illness develop, seek medical attention immediately.
Laceration: The Palodent® matrices are sharp and may cut unprotected skin. Handling matrices with gloved hands is recommended. In the event a clean matrix causes a skin laceration, employ usual first aid measures. Discard and properly dispose matrix.
2. Used matrices and accessories are contaminated. Follow proper exposure control plans. In the event of an accidental exposure, follow recommended post-exposure procedures. Following use of the Palodent® System, proper handling and disposal techniques are required. Used matrices are contaminated and sharp. When handling or disposing, use precautions as when handling or disposing, other contaminated sharps.
3. Care should be taken to have good control on any instrument to protect patients from injury.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Directions For Use. Any use of this product inconsistent with the Directions For Use is at the discretion and sole responsibility of the practitioner.
2. Care must be taken when placing the BiTine® ring, due to the spring nature of the ring. Ensure ring is fully seated in forceps (see Step-by-step Instructions).
3. Excessive pressure on matrix or retaining ring in an apical direction may result in injury to periodontal tissues.
4. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.

2.3 Adverse reactions

1. Product misuse may cause failure resulting in respiratory, digestive, eye, skin, or mucosa injury. (See Contraindications, Warnings and Precautions)

2. Failure to properly employ accepted handling and disposal techniques may present biohazard. (See Warnings and Precautions)

2.3 Storage conditions

Store at temperatures between 10-24°C/50-75°F, away from moisture or excessive humidity.

3. STEP-BY-STEP INSTRUCTIONS

1. Pre-operative sterilization

Matrices are intended for single use. If desired, these components may be steam autoclaved one time (134°C, 3 minutes minimum holding time, unwrapped, 12 minutes minimum holding time, wrapped) prior to clinical use. BiTine® rings and forceps are not supplied sterile, and should be steam autoclaved prior to first and each subsequent use (see Hygiene section).

2. Place rubber dam or other suitable isolation technique.

3. Pre-operative ring placement (optional)

- To effectively use a contoured sectional matrix, there must be enough interproximal space to position the matrix without force or distortion. The BiTine® ring, when placed prior to and left on during, cavity preparation, provides the necessary separation with no additional chair time.
- To place the round BiTine® ring, grasp it at its widest diameter with a rubber dam forceps, spread the forceps to expand ring only enough to clear tooth height of contour. Do not over-expand. Over-expansion may result in permanent ring deformation. The locking slide on the forceps may be engaged to facilitate carrying and placement of the ring.
- Place the ring with one tine in each interproximal space from buccal to lingual adjacent to the surface to be restored. The direction the ring points, mesial or distal, does not matter, though placing it so that the ring points mesially is usually more easily achieved due to easier access.
- To place the elongated BiTine® II ring, grasp it so that the rubber dam forceps engages the ring in the "u-shaped" depressions above the prongs. Place as described for the round ring. Either ring may be used to provide the necessary separation when applied before the preparation is initiated.
- With the tines properly positioned interproximally, release the locking slide, and slowly release the forceps tension. When ring is securely seated, release the forceps completely and withdraw forceps.

Small parts under tension

- Be sure ring is securely seated in the forceps before introducing into the patient field.
- Be sure ring is securely seated apically before withdrawing forceps.

4. Cavity preparation

Complete cavity preparation as usual. Once the preparation is complete, remove the ring.

5. Matrix selection

The standard Palodent® contoured sectional matrix is recommended for most applications. The Mini matrix is designed for use with adolescents and patients with poorly erupted posterior teeth. It is also ideal for the direct-indirect technique. The Plus matrix is designed for larger cavities occlusogingivally. It augments a contoured matrix with a patented adjustable apron gingivally and an adjustable continuation of the contour at the marginal ridge.

6. Application of matrix

- Once the tooth is prepared and the ring removed, finger roll the matrix to the approximate tooth circumference.
- Grasp the matrix with the notch toward the occlusal for the standard matrix, with the dot toward the occlusal for the Mini matrix and with the longer flap toward the apex for the Plus matrix, with a forceps or cotton pliers. If the ring was placed properly before tooth preparation, then the band should slide easily into place adjacent to the space to be restored from the occlusal approach.
- Close the gingival margin by placing an anatomic wedge.
- Replace the BiTine® ring in the proximal plane, engaging the matrix and wedge, for stability against tooth structure. The tines may be sitting on the wedge, in front of the wedge, or behind the wedge.
- Lightly burnish the matrix against the adjacent tooth both buccally and lingually to further adapt it to form proper contours. **NOTE:** Over burnishing may distort matrix and compromise contours.

Technique tip: When restoring an MOD preparation, two round rings may be placed facing opposite directions, one mesially and one distally or both pointing in the same direction, usually mesially, one over the other placing the mesial band and ring first. If placed pointing in the same direction, it is often easier to use one round BiTine® ring and one elongated BiTine® II ring. In this scenario, the round ring and band should be placed first mesially, then the elongated ring and band placed distally. Placing one ring with the tines pointing occlusally and one with the tines pointing gingivally also facilitates restoration of an MOD preparation.

7. Placement of restorative

Place restorative material and establish proximal contact as per usual technique. Light cure and/or allow material to reach initial set. Gross contouring may be accomplished with matrix in place.

8. Matrix removal

Remove ring(s), remove wedge(s) and remove sectional matrices. Finish contouring and verify appropriate proximal contact.

Technique Tips: It is advisable to "post-cure" light cured restorative materials following matrix removal from all surfaces, including proximo-facial and lingual line angles.

4. HYGIENE

4.1 Cleaning

The Palodent® system BiTine® rings may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.

4.2 Sterilization and maintenance

- Do not attempt to sterilize or reuse used sectional matrices. Properly dispose used matrices with other contaminated sharps.

Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

- The Palodent® System BiTine® rings should be steam autoclaved at 134°C for a minimum of 3 minutes. Rings may be autoclaved up to 700 times. Use of other sterilization procedures/methods is at the discretion and sole responsibility of the practitioner. It is not recommended to submerge Palodent® System components into disinfection or cold sterilization solutions.
- Sterilization of the matrices pre-use (if desired) may be accomplished (one time) by steam autoclaving, 134°C, 3 minutes minimum (unwrapped), 12 minutes minimum (wrapped).
- Allow BiTine® rings to dry prior to storage.

5. LOT NUMBER AND EXPECTED USEFUL LIFE

1. The following numbers should be quoted in all correspondence:

- Reorder Number
- Lot number on the package
- Copy of proof of purchase

2. The BiTine® rings may be used and reprocessed for use up to 700 cycles.

Manufactured by:
DENTSPLY Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
Tel.: 1-302-422-4511
www.dentsply.com

DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
Germany
Tel.: 49-7531-583-0
www.dentsply.de

Distributed by:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada

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