DIRECTIONS FOR USE – ENGLISH

For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

SureFil® High Density Posterior Restorative is a visible light activated, radiopaque composite restorative material designed for stress-bearing posterior restorations of primary and permanent teeth. The handling characteristics allow rapid, bulk placement and assist with establishment of properly located and shaped proximal contact areas.

1.1 Delivery forms

SureFil® Restorative material is available in:

- Individual dose, blister pack cups
- 3g syringes
- 3 shades: A-range, B-range, C-range

1.2 Composition

Dimethacrylate resin; TEGDMA dimethacrylate multifunctional polymethacrylate; Camphorquinone (CQ); Ethyl-4(dimethylamino)benzoate; Butylated hydroxy toluene (BHT); UV stabilizer; Silanated barium-boron-fluoro-alumino-silicate glass (mean particle size < 1 µm); Highly dispersed silicon dioxide (particle size 10-20 nm); Fluorescent agent; Iron oxide pigments; Titanium dioxide

1.3 Indications

SureFil® Restorative Material is indicated as a direct restorative for cavity Classes I, II in posterior teeth.

1.4 Contraindications

SureFil® Restorative material is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins.

1.5 Compatible adhesives

SureFil® Restorative material is used following application of a suitable dentin/enamel adhesive and is chemically compatible with conventional methacrylate-based dentin/enamel adhesives including DENTSPLY adhesives designed for use with visible light cured composite restoratives (see complete directions for use of selected adhesive).

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

SureFil® Restorative material contains polymerizable methacrylates which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

2.2 Precautions

- This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- The SureFil® Restorative material syringe should be tightly closed immediately after use.

Interactions:

- Eugenol- and hydrogen peroxide containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.

- SureFil® Restorative material is light-cured material. Therefore, it should be protected from ambient light. Proceed immediately once the material has been placed.
- If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.

2.3 Adverse reactions

Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous membranes: inflammation, edema, sloughing. See Warnings.

2.4 Storage

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Keep out of direct sunlight and store in a well ventilated place at temperatures between 10°-24°C/50°-75°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS

Direct restoration

1. **Shade Selection**
Before selecting the shade, teeth should be clean, hydrated and free of extrinsic material or stain. SureFil® Restorative Material shades, designated A, B and C, are slightly lighter than A2, B1 and C2 respectively. This deliberate slight mismatch permits visualization of the restoration margins during placement and finishing.
2. **Cavity Preparation**
 - Prepare the cavity so that no residual amalgam or restorative material is left.
 - Rinse surface with water spray and carefully dry it with air spray. Do not desiccate the tooth structure.
 - Use a dental dam or cotton rolls to isolate the cavity from contamination.
3. **Placement of Matrix**
The use of a sectional (e.g. Palodent®/Palodent® Plus Sectional Matrix System) or thin matrix band (e.g. AutoMatrix® Retainerless Matrix System) and subsequent burnishing of the matrix band will improve final interproximal contact and contour. Pre-wedging/ring placement is recommended.
4. **Pulp protection, Tooth Conditioning/Dentin Pretreatment, Adhesive Application**
Refer to adhesive manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of SureFil® Restorative material.
5. **Placement of SureFil® High Density Posterior Restorative**
 - Syringe**
 - Remove the cap.
 - Turn the handle of the syringe slowly in a clockwise direction and dispense the necessary amount of the material onto a mixing pad.
 - Point the front tip of the syringe upwards and turn the handle anti-clockwise to prevent oozing of the material.
 - Immediately close the syringe with the cap.
 - Protect the restorative material on the mixing pad against light.

Individual-dose Cups

Remove protective covering. Place material onto a clean mixing pad or dappen dish. Use the unit dose cup, inverted over remaining material, to protect SureFil® Restorative Material from premature polymerization via ambient light while current increment is adapted and light cured.

- 5.1 **Optional Flowable Liner:** Placement of a compatible flowable liner (available separately) prior to placement of SureFil® Restorative material is optional. Follow manufacturer's directions for use.
- 5.2 **Load nonstick, plastic insert type amalgam carrier or cut a suitable section with a sharp instrument.** Carry material to the preparation with the carrier or suitable placement instrument.
 - 5.2.1 **Class II:** Fill the proximal box in bulk (up to 5mm), level with the pulpal floor. Firmly adapt material to cavity walls and stabilize matrix band by stepping a clean non serrated condenser (plugger) from middle to edges, reducing incorporated air and eliminating voids. Use an amalgam condenser to pack the material. Light cure for 40 seconds (see Curing, Step 6). The remainder of the preparation is filled in bulk increments (up to 5mm).
 - 5.2.2 **Class I:** SureFil® Restorative Material may be placed and adapted in up to 5mm increments, light curing each increment. (See Curing, Step 6)
- 5.3 **Prior to light curing, contour and shape the final increment with the operator's choice of clean carving and burnishing instruments.** It is strongly recommended that the instrument used to contour cavosurface margins (occlusal and proximal) should be lubricated from time to time with a thin coat of residual adhesive to ensure optimal

marginal adaptation. **Important Technique Tip:** To minimize finishing time, contour marginal ridge first before occlusal anatomy. Carve margins and anatomy to final form. SureFil® Restorative Material resists slumping, allowing carving of the majority of anatomical form prior to visible light curing.

6. Curing

Light cure each area of the restoration surface with a suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470nm. Minimum light output must be at least 550mW/cm² exposure for at least 40 seconds. Some advanced performance curing units have been shown to cure 5mm increments of SureFil® Restorative material in 30 seconds. Refer to curing light manufacturer's recommendations for compatibility and curing recommendations. The SureFil® Restorative material should be additionally exposed to the curing unit through the proximal, lingual, and buccal enamel walls following matrix removal for the recommended time.

Inadequate polymerization due to insufficient curing

- Check compatibility of curing light.
- Check curing cycle.
- Check curing output before each procedure.

7. Finishing and Polishing

- Contour the restoration using finishing burs or diamond finishing instruments. Additional finishing is recommended by the use of Enhance® Finishing System. See manufacturer's complete directions for use.
- To achieve a high luster on SureFil® Restorative material, it is necessary to complete the polishing. PoGo® One Step Diamond Micro-Polisher System and/or Prisma®•Gloss™ Composite Polishing Pastes are recommended. See manufacturer's complete directions for use.

4. HYGIENE

Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

4.1 Cleaning and disinfection of the syringe

- To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination.
- Disinfect contaminated syringes with a water-based hospital-level disinfection solution according to national/local regulations.
- Repeated disinfection may damage label. **NOTE:** Vigorous wiping can destroy the label. Wipe syringe gently.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard uses: "YYYY/MM"
2. The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot number
 - Expiration date

©2012 DENTSPLY International. All Rights Reserved.

Manufactured by:
DENTSPLY Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
Tel.: 1-302-422-4511
www.dentsply.com

DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
Germany
Tel.: 49-7531-583-0
www.dentsply.de

Distributed by:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada

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