IRM® Caps™

DENTSPLY

Intermediate Restorative Material

DIRECTIONS FOR USE – ENGLISH

For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

IRM® Intermediate Restorative Material is a polymer-reinforced zinc oxide-eugenol composition restorative material designed for intermediate restorations intended to remain in place for no longer than one year. IRM® Intermediate Restorative Material may also be used as a base under restorative materials and cements that do not contain resin components.

1.1 Delivery forms

- Pre-dosed, encapsulated material for machine mixing, minimum deliverable volume 0.015mL

Alternatively available as:

- A powder/liquid version for hand mixing

1.2 Composition

Powder: Zinc oxide; Poly-Methyl Methacrylate (PMMA) powder; Zinc Acetate; Pigment

Liquid: Eugenol; Acetic acid

1.3 Indications

- Intermediate (temporary, less than 1 year) direct restorative filling material.
- Base material under non-resin-containing permanent restorations.

1.4 Contraindications

- Use with patients who have a known hypersensitivity or severe allergic reaction to eugenol, or any of the components.
- Direct application to dental pulp tissue (direct pulp capping).
- Base under resin-containing adhesives, restoratives or cements.

2. SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Liquid contains eugenol. Powder contains polymethacrylate. Components and mixed material may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible individuals.
 - **Avoid Eye Contact** to prevent irritation and possible corneal damage. In case of eye contact, rinse immediately with plenty of water and seek medical attention.
 - **Avoid Skin Contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - **Avoid contact** with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
2. Improved or reinforced Zinc oxide – eugenol products such as IRM® Material are not suitable for use as root canal filling materials as inconsistent, incomplete or delayed resorption may occur.

2.2 Precautions

- This product is intended to be used as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
- Insufficient data exist to support the use of IRM® material in a root canal as a cement, dressing, medicament or repair material.
- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Increased temperature and/or humidity reduce available work time and accelerate set time. Material will set faster in the mouth than at the bench. Do not use bench cure rates to judge mouth cure rates.

- Do not open capsule prior to trituration.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse.
- Use in a well-ventilated area.
- Slow speed amalgamators or slow speeds on amalgamators are not recommended.
- **Interactions**
 - Do not use this product in conjunction with (meth)acrylate based dental materials since it may interfere with the hardening of the (meth)acrylate based dental materials.

2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with Mucous Membranes: inflammation. (See Warnings)

2.4 Storage

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store in a well ventilated place at temperatures between 10-24°C/50-75°F.
- Allow material to reach room temperature prior to use.
- Do not store with or near materials containing (meth)acrylate resin components.
- Protect from moisture.
- Do not freeze.
- Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Cavity preparation

1. Isolate the field of operation in preferred manner. Rubber dam is recommended.
2. Cavity preparation must provide for the mechanical retention of the product.
3. If caries is present, completely excavate using low-speed and/or hand instrumentation.
4. Place matrix and wedge where indicated.
5. Gently dry preparation with cotton pellet. Avoid desiccation.

3.2 Pulp protection

In deep cavities cover the dentin close to the pulp (less than 1mm) with a hard setting calcium hydroxide liner (Dycal® Liner, see complete Directions for Use).

3.3 Mixing IRM® Caps™ material

1. To activate the IRM® Caps™ capsule, hold vertically, grasp bottom half and FIRMLY TIGHTEN the top (see Fig.1). As top is tightened, you will feel a "snap" as the liquid is released. Continue tightening until you are sure it cannot be seated further.
2. This complete tightening forces the liquid into the mixing chamber, thereby assuring proper handling and physical characteristics. Without delay, insert capsule (bottom first) into the amalgamator arms (see Fig. 2). Figure 4 gives approximate mixing times for various amalgamators. These are suggested ranges; however, your clinical experience in obtaining a good mix consistency with your unit remains the best guide to mixing time. Use the first five caps to establish the best mix time. (See 3.4 Guidelines for Optimal Mix)
3. Remove press cap (see Fig. 3) to gain access to the mixed IRM® Material.

Fig. 1


Fig. 2

Fig. 3

Fig. 4: Approximate Mixing Time (seconds)

AMALGAMATORS	SETTING	SECONDS
ProMix™	High Speed (Rabbit)	12
Vari-Mix®-III	High Speed	12
Silamat® Model C		8
Wig-L-Bug®*	High Speed	10
Wig-L-Bug™ Model 80	High Speed	10
Pelton Crane*		30
Capmaster*	Not Recommended	

*Not registered products of DENTSPLY International

3.4 Guidelines for optimal mix

1. Clinical experience with your amalgamator is the best guide to mixing the IRM® Caps™ capsule. Variation in amalgamators (even from the same manufacturer) will occur and the trituration speed of amalgamators will vary with the line voltage on which they are operated. High-speed amalgamators are recommended for best results.
2. If a dry, stiff or rubbery mix is obtained, reduce trituration time, keep speed the same.
3. If a sticky and soft mix is obtained, increase trituration time, keep speed the same. Sticky mixes should be allowed to set a bit before application, or rolled in IRM® powder to eliminate stickiness.
4. If trituration results in a non-coalesced mass (i.e., small balls), increase trituration time, keep speed the same.
5. If the mixed mass of IRM® material does not appear in the press cap (see Fig. 3) then either replace the press cap and mull for 5 seconds or insert a small ended plugger into the cap and remove the mixed mass.

3.5 Placement of mixed IRM® material

Mixed material may be carried to the preparation with suitable placement instrument. Adapt, contour and shape with appropriate instruments.

1. Use as Intermediate Restorative Material

- Trim excess material. Material may be smoothed by lightly burnishing. Initial set occurs approximately 5 minutes from start of mix. Increased temperature, humidity, and powder/liquid ratio may accelerate set.
- Upon set, remove matrix band if present. Check and adjust occlusion as needed. Additional carving or occlusal adjustment may be accomplished with a round bur. Avoid gouging material with sharp burs at high speed. Additional polishing is not recommended.
- Periodically evaluate the intermediate restoration until final restoration is indicated. The intermediate restoration should be removed and final restoration completed within 1 year of placement.

2. Use as Protective Base

- This product is a protective base for restorative material not containing polymeric resin components (e.g. final restoration of amalgam, gold).
- **NOTICE: Possible interactions** – The product should not be used in conjunction with resin-based adhesives, varnishes or restoratives since eugenol may interfere with hardening and cause softening of (meth)acrylate based dental materials.
- Apply mixed material to desired dentin surfaces and depth.
- Following completion of set (approximately 5 minutes from start of mix) cavity refinement may be completed.
- Complete restoration per restorative material manufacturer's directions for use.

4. HYGIENE

Cross-contamination

- Do not reuse single use products.
- Dispose in accordance with local regulations.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard is used: "YYYY/MM."
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date

Manufactured by
DENTSPLY Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
Tel: 1-302-422-4511
www.dentsply.com

Distributed by
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada

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DENTSPLY DeTrey GmbH