Calibra® Veneer

Esthetic Resin Cement

DENTSPLY

DIRECTIONS FOR USE – ENGLISH

For dental use only. USA: Rx only.

1. PRODUCT DESCRIPTION

Calibra® Veneer Esthetic Resin Cement is a visible light cured, high strength resin cement compatible with numerous dentin/enamel adhesive systems. It is to be used following the application of a suitable dentin/enamel adhesive to adhesively bond and lute indirect restorations to tooth structure.

1.1 Delivery forms (Some delivery forms may not be available in all countries)

Calibra® Veneer Cement and Try-In pastes are available in:

- small syringes for easy dispensing
- 5 shades: translucent, bleach, light, medium, and opaque

1.2 Composition

- Calibra® Veneer Esthetic Resin Cement Base: Dimethacrylate Resins; Camphorquinone (CQ) Photoinitiator; Stabilizers; Glass Fillers; Fumed silica; Titanium Dioxide; Pigments. Particles of inorganic filler range from 0.02 to 1.3µm, total filler 44.9% by volume.
- Calibra® Veneer Esthetic Resin Cement Try-In Paste: Glycerine; Fumed Silica; Titanium Dioxide; Pigments
- Calibra® Silane Coupling Agent: Acetone; Ethyl Alcohol; Organo Silane

1.3 Indications

Adhesive cementation of ceramic, porcelain, composite veneers.

1.4 Contraindications

Calibra® Veneer Cement is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins or any of the components.

1.5 Compatible adhesives

Calibra® Veneer Cement is used following application of a suitable dentin/enamel adhesive and is chemically compatible with conventional methacrylate-based dentin/enamel adhesives including DENTSPLY adhesives designed for use with visible light cured composite restoratives (see complete Directions for Use of selected adhesive).

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other sections of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Calibra® Veneer Cement contains polymerizable methacrylate monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
 - Avoid Eye Contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - Avoid Skin Contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and alcohol and wash well thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - Avoid Contact with Oral Soft Tissues/Mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
2. Calibra® Silane Coupling Agent contains acetone. Do not breathe vapor.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Insufficient data exist to support the use of Calibra® Veneer Esthetic Resin Cement as a core build up, filling material, base or liner.
3. Insufficient data exist to support the use of Xeno® Self-Etch Adhesives (all versions) with Calibra® Veneer Esthetic Resin Cement for veneer cementation light cured.
4. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
5. Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
6. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
7. The syringe cannot be reprocessed. To prevent the syringe and from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringe is handled with

- clean/disinfected gloves. As additional precautionary measure, syringe may be protected from gross contamination but not from all contamination by applying a protective barrier.
8. Containers should be tightly closed immediately after use. Syringe should be tightly closed by replacing the original cap immediately after use.
9. Calibra® Veneer Cement should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause syringe rupture.
10. Use only in well ventilated areas.
11. Flammable: Calibra® Silane Coupling Agent contains acetone. Keep away from sources of ignition.
12. Interactions:
 - Eugenol containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
 - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.
 - The Calibra® Veneer Cement is light cured material. Therefore, it should be protected from ambient light. Proceed immediately once the material has been placed.

2.3 Adverse reactions

Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous membranes: inflammation (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Keep out of direct sunlight and store in a well ventilated place at temperatures between 2°-24°C/35°-75°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS FOR USE

- I. Following removal of the temporary restoration and any remaining temporary cement, clean enamel and dentin as directed by the adhesive manufacturer's directions.
- II. Check the fit and esthetics of the restoration. **Technique Tip:** Occlusal adjustment of veneers is best accomplished after final cementation.

3.1 Try-in Paste (optional)

- Dispense the appropriate shade of try-in paste from the syringe onto a clean mixing pad or glass slab. Load paste onto internal surfaces of the restoration and gently seat onto preparation. Clean excess with a cotton pellet and/or blunt explorer. Shades may be blended to achieve optimum esthetics. **Important Technique Tip:** Factors such as restoration shape, texture, thickness and shade(s) of the underlying and surrounding tooth structure contribute to final perceived shade. Try-in paste is a guide for cement shade range selection only. **NOTE:** The try-in paste will not polymerize, thus offers unlimited work time.
- Once fit and esthetics are verified, thoroughly rinse try-in paste from restoration and preparation surfaces using water.

3.2 Treatment of the restoration: Ceramic/Composite Restorations

Follow the dental laboratory or restoration manufacturer's instructions for pre-treatment, if required. Restorations designed to be silanated or if the internal silanated surface has been disturbed during try-in, apply Calibra® Silane Coupling Agent as outlined below:

- Prior to applying the silane, clean the internal surface of the restoration with Caulk® 34% Tooth Conditioner Gel (34% phosphoric acid) (see complete Directions for Use). Apply for 30 seconds and rinse thoroughly with water for 20 seconds. Air-dry the internal surface of the restoration with clean, dry air from a dental syringe.

Danger of injury due to excessive force

- Apply slow and steady pressure on the syringe plunger
 - Do not use excessive force – syringe rupture may result.
- Attach supplied needle tip to end of the Calibra® Silane Coupler syringe. Gently apply pressure to syringe plunger. DO NOT USE EXCESSIVE FORCE. Silane agent should express easily, one drop at a time. Directly apply the silane to the clean internal surface of the restoration. Air-dry the silanated surface of the restoration with clean, dry air from a dental syringe or allow to air dry. **NOTE:** Avoid application of the silane to the external surface of the restoration. Contact with the external surface will cause bonding of the cement to the external surface, complicating clean-up and necessitating removal of the glazed portion of the external surface of the restoration.
 - If silanated surface becomes contaminated, clean surface with alcohol or acetone, allow to air dry and repeat application as outlined above.

3.3 Tooth conditioning/dentin pretreatment/adhesive application

- Following restoration try-in, rinse preparation thoroughly with water spray and air dry.
- Proceed with adhesive application as directed in the complete directions for use.

Delayed or reduced adhesion due to adhesive incompatibility

- Select compatible light cure adhesive
 - Follow selected adhesive directions for use
- Always refer to and follow adhesive manufacturer's complete directions for use for selection and application of the appropriate light cure adhesive. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately following adhesive application to cementation with Calibra® Veneer Cement, light cured base as indicated in the appropriate section below.

3.4 Veneer cementation – light cured

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment/Adhesive Application.

- Adhesives compatible with light cured, CQ initiated methacrylates are recommended.
- If veneer restorations exceed 1.0mm thickness, are heavily opaqued, or otherwise impede the transmission of light, operator is urged to select Calibra® Esthetic Resin Cement and follow the Dual Cured cementation technique, mixing base with catalyst (available separately).
- Apply adhesive to the internal bonding surface of the restoration if directed by adhesive manufacturer's Directions for Use.

Cementation technique

- Dispense the desired shade of Calibra® Veneer Esthetic Resin Cement base paste from the syringe directly onto the veneer. Protect cement from exposure to light.
- Seat the loaded veneer in place. Remove excess from the gingival margin with a blunt instrument. Briefly (10 seconds or less) light cure the gingival portion only to tack restoration in place. Remove any excess from proximal and lingual margins. **Technique Tip:** Placement of mylar strips between preparation and adjacent teeth prior to seating veneer aids in isolation and in excess cement clean up. After removing gingival excess and tacking veneer in place, remove interproximal excess cement by pulling mylar strip towards the facial, i.e., from tooth to restoration surface.

3.5 Curing

Inadequate polymerization due to insufficient curing

- Check compatibility of curing light.
 - Check curing cycle.
 - Check curing output before each procedure.
- Using a visible light* curing unit designed to cure CQ initiated methacrylates (spectral output including 470nm), with a minimum output of 550mW/cm², light cure 20 seconds each from the buccal, lingual, and interproximal aspects.
- Following the light cure check and adjust occlusion as necessary. Proceed to Finishing and Polishing section.
- *Refer to curing light manufacturer's recommendations for compatibility and curing recommendations.

3.6 Finishing and polishing

- Removal of resin flash is best accomplished with the Enhance® Finishing System of points, cups and discs. The Enhance® System will remove flash and finish restoration margins without removal of or trauma to the enamel.
- Polish final restoration using Enhance® PoGo® Polishing system or Prisma® Gloss™ Polishing Paste and Prisma® Gloss™ Extra Fine Polishing Paste. (See complete Directions for Use supplied with polishing product chosen)

All shades of Calibra® Veneer Cement are radio-opaque, with a 1mm radiopacity equivalent to 1.77mm radiopacity of aluminum. Aluminum has a radiopacity equivalent to that of dentin. Thus 1mm of material having a radiopacity equivalent to 1mm of aluminum has a radiopacity equivalent to that of dentin.

4. HYGIENE

Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.

4.1 Cleaning and disinfection

- To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid syringe contamination. The use of protective barriers is an additional precautionary measure against gross contamination but not against all contamination.
- Handling with clean/disinfected gloves, syringes may be disinfected with a water-based hospital-level disinfection solution according to national/local regulations.
- Repeated disinfection may damage label. **NOTE:** Destruction of the label by vigorous wiping. Wipe syringe gently.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard is used: "YYYY/MM"
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date

Manufactured by
DENTSPLY Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
Tel: 1-302-422-4511
www.dentsply.com

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DENTSPLY DeTrey GmbH