DENTSPLY
TRUBYTE

Lucitone 199[®] Repair Material

Résine de Réparation
Materiale per riparazioni
Reparaturmaterial
Material para reparaciones

DIRECTIONS FOR USE
CONDITIONS D'EMPLOI
ISTRUZIONI PER L'USO
VERARBEITUNGSANLEITUNG
INSTRUCCIONES DE EMPLEO

DENTSPLY
TRUBYTE

Manufactured by:
DENTSPLY International Inc.
570 West College Avenue
York, PA 17405-0872
1-800-786-0085

www.trubyte.dentsply.com or www.dentsply.com

EC REP

DeguDent GmbH
Postfach 1364 • 63403 Hanau
GERMANY
Telefon +49/6181/5950

580895 Rev. F (06/09)

R_xonly

CE

Indications for Use: LUCITONE 199® Repair Material is indicated primarily for repair and relining of removable dentures. This material may also be used for other applications in the dental laboratory where a chemical cure acrylic resin may be used, such as baseplates. The powder is available in all four Lucitone 199 shades and is formulated for use only with DENTSPLY Repair Material Liquid.

Contraindications: Lucitone 199 Repair Material is contraindicated for patients and users with a history of allergic reaction to methyl methacrylate monomer.

Warnings:

1. Lucitone 199 Repair Material contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical advice.
2. Avoid inhalation or ingestion. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Liquid contact with eyes may cause possible corneal damage. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.

Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.

Ingestion: Contact your regional Poison Control Center immediately.

Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician. Wash skin with soap and water.

Precautions:

1. Users with special skin problems, cuts, or abrasions should wear protective gloves.
2. When grinding prosthodontic resins, proper ventilation, masks and vacuum systems should be used.
3. Store at 60°-80°F, away from moisture and direct sunlight. The liquid contains unsaturated monomers which may polymerize prematurely if stored at excessively high temperatures and/or in direct sunlight. The powder contains an organic polymerization initiator which may degrade if stored at excessively high temperatures.


4. DENTSPLY® Repair Material Liquid and Lucitone 199 Repair Material uncured resin are hazardous materials. Dispose of in accordance with Federal, State, and Local regulations.
5. DENTSPLY Repair Material Liquid contains methyl methacrylate monomer, a flammable liquid with a flash point of 10°C (50°F). Keep away from heat, sparks, and open flame.
6. Use liquid in a well ventilated area. Replace cap when not in use.
7. Do not pack denture resin after work time has elapsed or material has become rubbery.
8. When using either the dusting or mixed-dough method, L199 Repair Material must be cured in a pressure curing unit.

Adverse Reactions:

1. Corneal damage, headache, nausea, and vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precaution statements.)
2. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthetics in warm water for several days.
3. Particulates will be generated acrylic resins. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.


STEP BY STEP INSTRUCTIONS:

A. REPAIRING INCOMPLETE FRACTURES

Small crack is running about halfway into the palate from posterior border.

B. REPAIRING COMPLETE FRACTURES

The complete fracture is assembled and held in position with splints and sticky wax; cast is poured. For a denture having severe undercuts, a cast may be poured using Dentsulate™ Insulating Putty.

C. ADDING TOOTH AND CLASP

Denture is seated on tinfoil or AL-COTE® Separator coated cast. The clasp is in place and the tooth is resting in the plaster matrix.

D. POST DAMMING

One-half inch of the posterior border of the denture is roughened. The extension is made with LUCITONE 199 Repair Material.

E. REPLACING BROKEN OR DISLODGED TEETH

Upper right central and lateral are fractured.

F. IMMEDIATE PARTIAL (TEMPORARY)

Tooth to be extracted is cut from the cast. Design of the partial is outlined on the cast.

PREPARATION OF REPAIR AREAS

The preparation of repair area may vary slightly based on the circumstances encountered. The following recommendations are general in application. When repairing a broken denture, hold the fractured parts in position with splints and sticky wax while pouring the stone cast. Next, apply AL-COTE® Separator to the cast, open the repair area to approximately ¼ inch, and bevel the edges of the opening with a narrower width on the tissue side of the repair.

Complicated dovetailing weakens the repair and wastes time.

Edges of repair areas should be smoothed with pumice or with a sharp knife to avoid a frosty demarcation line.

LUCITONE 199 Repair Material may be used with either the dusting or packing (as for relining) method.

DUSTING PROCEDURE

When using the dusting method, LUCITONE 199 Repair Material must be cured in a pressure curing unit.

Dust the powder in successive thin layers and moisten each layer with DENTSPLY Repair Material Liquid. The repair area should be slightly over-contoured for best results.

Porosity and/or surface blanching may occur if a pressure curing unit is not used.

Submerge the required appliance (and cast as appropriate) in a warm water bath (110°F) under a minimum of 15 psi pressure for 10 minutes. The result is maximum density in the repaired area.

PACKING PROCEDURE

It may be desirable to invest (flask) and pack LUCITONE 199 Repair Material for procedures such as post damming, fabricating temporary partials and processing indirect relines.

When packing, follow the mixing and curing methods described for indirect relining procedures.

INDIRECT RELINE PROCEDURES

The impression is obtained using an old denture as a tray. The periphery is trimmed and large undercuts eliminated before taking the impression.

LUCITONE 199 Repair Material is highly recommended for restoring a comfortable fit to poorly fitting dentures by the indirect method.

Take an impression in the relieved denture. Then mount the stone cast with the denture in position on an articulator or duplicator to preserve the tooth/cast dimension.

Remove the denture and clean out most of the impression material, leaving sufficient impression material to serve as automatic rests. This will help assure a positive reseating on the cast. Alternatively, use matrices for the tooth/cast relationship.

If indicated, scrape or sand the high spots and periphery to further relieve the denture. Smooth all surfaces. Reseat the denture on the stone cast. Apply wax as needed to areas that have been ground away.

Invest the case in the customary manner. Heat flask in hot water, separate and remove remaining impression material and flush out the wax. Apply AL-COTE Separator to the exposed investment and on the cast. Cool the flask to room temperature.

Make a mix of LUCITONE 199 Repair Material using about 15 cc of powder and 5 ml of liquid - a ratio of 25 grams of powder to 17 ml liquid. **In about 14-20 minutes, the mix will reach a good packing consistency, with a softer texture than that of an ordinary denture resin.**

Wet the tissue side of the denture with liquid and pack the mix of Lucitone 199 Repair Material. No trial packing is necessary. Maintain pressure in clamp or press for 15 minutes, then deflask the case. It is not necessary to use water or a pressure pot with this packing method.

Completed denture after relining with LUCITONE 199 Repair Material by indirect method.

FINISHING AND POLISHING

Grind off all excess repair material (down to the denture surface level) with an appropriate trimmer or bur. Chisel away flash around the necks of the teeth when repairing teeth.

To polish, use brush and rag wheels with pumice, followed by a denture high shine.

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