Luci-Sof®

Durable Soft Denture Liner System

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DIRECTIONS FOR USE

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Indications for use:

The Luci-Sof Soft Denture Liner System is intended for use as a long-term resilient lining material for removable dentures.

The Luci-Sof Denture Liner System consists of two components:

1. Luci-Sof Soft Denture Liner: a heat polymerizable, one component silicone putty that is based on dimethylpolysiloxane, dimethyl, methyl hydrogen siloxane copolymer; reinforced with amorphous silica. Benzoyl peroxide is the initiator for polymerization.
2. Luci-Sof Bonder: a heat polymerizable methacrylate monomer liquid and methacrylate polymer blend.

Contraindications:

1. Luci-Sof Bonder is contraindicated for patients and users with a history of allergic reactions to methyl methacrylate monomer.

Warnings:

1. Luci-Sof Bonder contains polymerizable monomers which may cause skin sensitization (allergic contact dermatitis) or other allergic reactions in susceptible individuals. Wash thoroughly with soap and water after contact. If skin sensitization or other allergic reactions occur, discontinue use. If dermatitis or other symptoms persist, seek medical advice.
2. Avoid inhalation or ingestion of Luci-Sof Bonder. High vapor con-

centration can induce headache, nausea, and irritation of eyes and respiratory system. Liquid contact with eyes may cause possible corneal damage. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.

Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.

Ingestion: Contact your regional Poison Control Center immediately.

Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician. Wash skin with soap and water.

Precautions:

1. Do not use Luci-Sof Soft Denture Liner on dentures that will result in denture base acrylic that is less than 2 mm in thickness.
2. Do not use the alternative long cure cycle (9 hours at 163° F) when lining a new denture.
3. Users with special skin problems, cuts or abrasions should wear protective gloves.
4. When grinding prosthetic resins, proper ventilation, masks and vacuum systems should be used.
5. Refrigeration of Luci-Sof Soft Denture Liner is recommended. Always bring to room temperature before use. Do not store next to amine containing chemicals (ie. acrylic repair liquid, adhesive

accelerator, etc.) as this may cause discoloration of the Luci-Sof material.

6. Store Luci-Sof Bonder below 80° F (27° C) and away from direct sunlight. May be refrigerated to extend shelf life.
7. Eugenol containing dental materials should not be used in conjunction with this product since they may interfere with hardening and bonding.
8. Unused Luci-Sof Soft Denture Liner should be polymerized prior to disposal.
9. Luci-Sof Bonder is a hazardous material. Dispose of in accordance with Federal, State, and Local regulations.

Adverse Reactions:

1. Corneal damage, headache, nausea, and vomiting may occur with exposure to methyl methacrylate monomer. (See Warnings and Precautions statements.)
2. Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.
3. Particulates will be generated when grinding prosthodontic materials. Eye, skin and respiratory irritation may occur if appropriate engineering controls are not used.

STEP-BY-STEP INSTRUCTIONS:

Relining An Existing Denture

Luci-Sof Soft Denture Liner System is a long-term silicone heat-cured soft denture liner material and bonder.

1. Pour a stone cast of the impression in the denture to achieve the desired contour.
2. Trim excess impression material, smooth with wax and finish to final shape around the periphery. **This preparation will help to minimize finishing later.**
3. Invest.
4. Open flask and remove wax with boiling water.
5. Remove impression material.
6. Wash with detergent and water and rinse with clean boiling water.
7. Grind the denture surface with a clean carbide cutter to remove tray adhesive and sufficient acrylic. **A clean surface is necessary to achieve a superior bond to the denture base.** Avoid knife edges. Grind flange border to a 90° butt joint.
8. Shake the Luci-Sof Reline Bonder well and generously coat the areas

of the denture to be bonded.

9. Wash brush and dampen brush with monomer and dry with a paper towel.
10. Allow the Bonder to dry for 15 minutes.
11. Trial pack the Luci-Sof Soft Denture Liner. Add material to borders and heel areas to achieve uniform flash.

12. Immediately process in a compressor for 2½ hours in boiling water. (Alternative cure method: 9 hours at 163° F followed by a 30 minute boil.) Cool 15 minutes on the bench and then place in room temperature water for 15 minutes.
13. Remove excess material with a scalpel and contour with a Faskut Luci-Sof Bur. (Large #904765, Medium #904764, Small #904763)
14. Polish with Tripoli and felt wheel at high speed.

Lining A New Denture

1. Flask waxed-up denture using conventional methods.
2. Boil out and wash with detergent and water, then rinse with clean boiling water.
3. Prepare a 2-3 mm shim on the cast.
4. Pack denture base material and close flask. Bench set denture base under pressure until a firm

5. consistency is achieved (60 min. for Lucitone 199® Denture Base). Bench set time may vary for other brands of denture base.
6. Remove shim and trial pack the Luci-Sof Soft Denture Liner. **Do not use Luci-Sof Bonder.** Underpack the first trial and then add material to heel areas and borders.
7. Immediately process in a compressor for 1½ hours at 163° F, followed by 2½ hours in boiling water. Cool 15 minutes on the bench and then place in room temperature water for 15 minutes.
8. Remove excess material with a scalpel and contour with a Faskut Luci-Sof Bur. (Large #904765, Medium #904764, Small #904763)
8. Polish with Tripoli and felt wheel at high speed.

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DENTSPLY International Inc.
DeguDent GmbH