



DENTSPLY International, Inc. is a multinational health care corporation with operations around the globe and a rich history of product excellence. For over a century, DENTSPLY International has been committed to providing the dental community with innovative, high-quality, cost-effective dental products, services and equipment. From humble beginnings in 1899, the company has grown to become the largest professional dental products manufacturer in the world.

DENTSPLY Ceramco® is one of the world's premier brands for crown and bridge restorative products and materials. Ceramco porcelains are renowned for their quality, natural appearance and durability, characteristics that give dentists confidence that their restorations will deliver the health and esthetic benefits their patients desire.

DENTSPLY International, Inc.
570 West College Avenue
York, Pennsylvania 17405-0872
1.800.243.1942
prosthetics.dentsply.com
dentsply.com



1 Vita® and 3D-Master® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co.

When everything comes together.

Beauty. Integration. Consistency. Ceramco® iC.



When everything comes together.

See what Ceramco® iC can do for your patients.

Exceptional Beauty.

Designed to be highly esthetic, the Ceramco® iC ceramic system allows the creation of exquisite porcelain restorations. With a full spectrum of 46 shades, including 16 A-D shades, 26 CC Shade Series shades, keyed to the Vita® 3D-Master®¹ Shade Series, and 4 bleach shades, almost any prescription can be filled for exceptional patient satisfaction.

Unique Integration.

Ceramco iC is the unique integrated ceramic system that offers a single set of porcelains, pressing materials and stains for single or multi-unit porcelain-fused-to-metal and press-to-metal restorations, as well as for all-ceramics including crowns, veneers, inlays and onlays. Therefore providing the ultimate in shade consistency and ensuring a reliable, life-like shade match on combination cases involving more than one type of substructure.

Shade Consistency.

The Ceramco iC system imparts consistent fluorescence throughout all shades. The colors and graded opacities of the system components are maintained within controlled ranges and with the vitality of the Ceramco iC opalescence, the optical characteristics of natural dentition can be duplicated consistently and repeatedly.



Consistent
Fluorescence



Vibrant
Opalescence

Features

- **Increased translucency** for enhanced beauty crowns
- **Integration** across all-ceramic, porcelain-fused-to-metal and press-to-metal restorations
- **Shade consistency** for combination cases
- **Improved flexural strength** of 90 MPa for PFM and 135 MPa for PTM and All-Ceramic, allowing for thin veneers to 0.3mm.

Benefits

- Exceptional esthetics and patient satisfaction
- Expanded range of indications
- Reliable shade results and high-quality restorations



RESTORATIONS CREATED AND PHOTOGRAPHY PRODUCED BY
J. GEBHARDT, CDT, VANGELES, L.A., CALIFORNIA.